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Maxillary central incisor restoration after bone reconstruction 5 year Results

A 43 year-old patient consulted after she had lost her central incisor (11) by accident. During the diagnosis, we observe that the vestibular bone table has resorbed but the height is preserved. We propose a treatment plan including a subepithelium-connective tissue graft followed, 3 months after, by an autologous bone grafting. The bone bloc is harvested from the retromolar region.

6 months later, a Ø 4 mm x L. 10 mm twinkon® implant is placed in one surgical time. After a healing period of 4 months, the final prosthesis is fabricated.

At 5 year follow-up, we observe a perfect stability of the soft tissues as well as the formation of a cervical cortical seems to show on X-rays.

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

Implant supported restoration of the molar regions of the upper and lower jaws

4 year results


A female patient, 59 year old, executive in a private company, consults because she is edentulous on the maxillary premolar sectors and mandible molar sectors. The patient suffers a non-insulino dependent diabetes. Furthermore, she cannot accept the idea of a removable overdenture and would like a fixed restoration alternative. After diagnosis, we propose her an implant supported prosthetic solution.

Necessary extractions are performed in sectors 2 and 3. The sockets are filled in the same session with BIOBank cancellous bone powder. A diagnosis, driven with radiologic and surgical guides, shows a


moderate atrophy of the sub-sinus region as well as the horizontal ramus. The patient accepts to receive an allogeneic bone graft which is performed in March 2011 by Dr. V. Duguet. On September 2011, short In-Kone® implants are placed in one session and buried during a 3 month osseointegration period. The healing screws are placed on a second surgical time and, on March 2012, the final single units prosthesis are fabricated. The restoration is fully completed on April 2012.

On August 2016, 4 years and 4 months later, the patient who is 65 year old is now retired and fully satisfied with the rehabilitation.

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


Multiple maxillary and mandibular agenesis - 2 year results

A 16-year old female patient with multiple agenesis ectodermal dysplasias, particularly at the maxillary comes to our office. She has already been surgically treated to correct a class III malocclusion.

The orthodontic translation of 23 has failed. Our diagnosis starts from this situation. We decide to wait until the patient has grown up to undertake widening parietal grafts on the posterior area before placing 6 EVL® implants in the maxillary and 2 in the mandible.

After healing, the three anterior teeth are removed, and 2 other implants are placed. A complete provisional bridge is fabricated and loaded onto the 6 posterior implants.

The final prosthesis is fabricated, screw retained on the anterior teeth and cemented on the posterior ones.

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Management of an implant failure in the aesthetic area

A 48 year old patient consults for the emergence of two implants in positions 21 and 22 are not compatible with the fabrication of implant supported crowns. Furthermore, the patient has a bridge on 11-12-13 and a crown on 23. In order to obtain a better aesthetic result, we decide to place two new implants in positions 12 and 22 together with a bridge to restore the 4 incisors.

First of all, we remove the inappropriate implants and realize a provisional bridge together with a connective tissue graft in order to shape the soft tissues contours.

Then, 8 weeks after healing, we place two In-Kone® implants in 12 and 22 and proceed to a guided bone regeneration (drill bone, bio-oss, Kreos membrane). After 4 months of osseointegration, we remove the 11 and improve the gingival profile thanks to a new connective tissue graft together with a second provisional bridge.


Finally, 8 weeks later, we fabricate the final prosthesis: a 4-element bridge cemented on abutments and two ceramic crowns on 13 and 23.

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


Multiple prosthetic restorations on implants 16 year results

This patient suffered agenesia of the teeth: 15, 14, 13, 23, 24, 25, 34, 35, 44 and 45. An orthodontic treatment kept the inter-dental spaces stable until the patient is 18. The milkteeth are removed and two small metal based overdentures are placed temporarily. The patient hardly stands these prostheses. In 2000, he decides to undertake an implant treatment. A sinus lift is needed in the area 20. A first surgery is performed under neuroleptic to place 7 EVL® implants as well

as the sinus lift. The implants are inserted with a minimally invasive protocol. The area 20 will receive implants 6 months later. The loading with the final prosthesis is achieved 6 months after implants placement. Composite occlusal wedges are fabricated on 36 and 46 in order to improve the bite.

16 years later the case is stable, without bone absorption. The prosthesis is still functional.

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


Single Unit restoration in aesthetic zone 5 year results

A 53 year old patient consulted because of a fistula on 12. After X-ray analysis, we decided to remove the tooth and to perform a subepithelial connective tissue graft in order to enhance the healing of the extracted site and to increase the width of the mucosa.


After a healing period of 3 months, an autogenous graft was performed. The bone was harvested from the retro molar region in accordance with Pr. Khoury's biologic grafting technique. An In-Kone® implant

has been placed 4 months after the bone grafting. The good primary stability allows us to immediately load the implant with a provisional restoration. The final restoration is fabricated after a 3 month osseointegration period.

The X-ray control after 5 years shows the bone remodelling around the neck of the implant and the stability of the peri-implant tissues.

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Full-arch upper jaw and posterior mandibular restorations

A male 61 years old consults us because his upper bridge is instable. In the meantime, he would like a complete rehabilitation in order to replace his removable overdenture at the lower jaw as well as to recover an acceptable aesthetic aspect. The diagnosis shows a terminal periodontal disease at the terminal stage on the upper jaw as well as on the molars of the lower jaw. The different prosthesis previously fabricated lead to the inversion of the occlusion curves together with an unsightly smile.

After the teeth extraction, the prosthetic project is validated thanks to a removable prosthesis. In spite of the curves correction together with the aesthetic improvement, the edges of the maxillary incisors are still

invisible under the upper lip at rest. We thus decide to lengthen further the teeth. A new trial is carried out to validate the definitive project.

The final prosthesis is fabricated on a CAD/CAM bar which is screw retained on 6 In-Kone® implants at the upper jaw. Two screw retained ceramic bridges are fixed on the posterior sectors of the lower jaw. Both In-Kone® 6 mm short implants and twinkon® 4 mm ultra-short implants are placed in the latter case because of the reduced bone height available. The emergence of the dental nerve on the left prevents us from placing an implant in 34. Therefore a bridge is proposed as an alternative.

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
Replacement of a central incisor 3 year results


This 28 year old patient consults in emergency after being hit during a rugby match. The clinical and X-ray diagnosis reveal a longitudinal root fracture of the 11. This tooth is 'condemned' and needs to be removed. The analysis of the different factors leads our therapeutic choice towards extraction and immediate implant placement and loading with a provisional crown. The protocol consists in the placement of a Ø 4mm In-Kone® implant in the palatal wall of the socket straight after tooth extraction. The vestibular gap is filled with bovine

bone graft material. The marginal gingiva is thickened with a connective tissue graft harvested in the maxillary tuberosity region. A screw retained provisional restoration is fabricated on a H3 mm temporary abutment. Four and a half months later, once the implant is osseointegrated, the impression is sent to the technician. A ceramic crown on a zirconia basis is fabricated and glued onto a titanium basis. Finally, the crown is tightened on the implant and a suitable aesthetic integration is obtained after 3 years.

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
Full-arch implant supported prosthesis on electro-galvanic copings

A female patient with a full edentulous maxillary and partial edentulous mandible is asking for a fixed rehabilitation on implants.

After clinical and X-ray diagnosis we propose the following treatment plan:

- Extraction of the residual mandibular teeth and immediate implantation with 6 In-Kone® implants.
- Immediate placement of a provisional prosthesis Brånemark type.
- After a healing period of 3 months, fabrication of a removable Preteau Zirconia prosthesis stabilized by electro-galvanic copings. Then, fabrication of a new removable overdenture on the maxillary in order to restore the occlusion as well as the aesthetic aspect.

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


Restoration of a central incisor 8 year results

A 22 year old female patient consults in August 2007. She suffers from agenesis on 22 which is corrected by an unsightly cemented bridge. The patient is in good health and would like to restore the missing tooth with an implant. The clinical diagnosis reveals a sufficient mesio-distal space but convergent crowns. This diagnosis is confirmed by X-rays. The vestibular palpation indicates that the bone volume is suitable with the placement of an implant. We address our patient to an orthodontist in order to slightly open the available space and straighten the adjacent teeth.

The patient who is using aligners comes back in June 2008. A diagnosis wax up is prepared and a surgical guide is fabricated. An EVL[®]S 3.3/14

mm is placed in November 2008 in one surgical time. The roll flap technique is performed in order to augment the vestibular soft tissues and to modify the biotype. A flat cover screw (0.4 mm) is placed to close the implant. 8 weeks post-op, a X-ray is performed and the healing screw is switched with a higher one (3 mm). The first impression takings are made in March 2009 in order to fabricate a temporary crown on a FM Tempo abutment. This temporary crown will be modified on 4 different steps in order to guide the soft tissues to the desired architecture. After a new impression taking of the implant and the remodeled soft tissues, a trans-screwed ceramo-metal crown is placed on October 2009. Clinical and X-ray controls are made every year.

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Restoration of two central incisors 3 year results

A 39 year old female patient consulted because she is not satisfied with her smile. After the diagnosis, we suggest to remove the already treated central incisors and to immediately place two In-Kone[®] implants loaded with provisional crowns. The latter are replaced by second provisional restorations two weeks post-op in order to guide the maturation of the peri-implant soft tissues to the desired architecture. 3 months later, the final restoration is delivered.

But, rapidly the patient is not satisfied with the color of the crowns when the sun reflects on them. After we control under polarized light, we can clearly identify the color distortion. Accordingly, we fabricate new crowns.


X-rays at three years show stabilized peri-implant tissues. We regret though a small scar which remains and which is probably due to the initial suturing. Today, the patient is fully satisfied.

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

Complete upper jaw restoration with immediate loading 7 year results

A 79 year old female patient consulted because of her bridge mobility. After diagnosis, we propose to our patient a treatment plan including extraction and immediate placement and loading of 8 In-Kone® implants within the same session.

Our operating procedure is based on a very precise surgical protocol built on a symmetric positioning of the implants in native bone.

From an aesthetic point of view, our protocol also supports the preservation of the patient mucosa together with the use of grafting combined with PRF. Long term repeatable results can be obtained under the condition that the protocol is strictly respected.


At 7-year post surgery, note the quality of the results from the aesthetic point of view as well as the health of the peri-implant tissues.

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



Lateral incisor restoration 4 month result

A 45 year old female patient is referred with a fractured hopeless lateral incisor #22 for implant therapy. The root is atraumatically extracted and a narrow 3.0 Implant x 13mm implant is then placed in the extraction alveolar in the proper 3D position. A tulip shape transgingival component with a Ø 3,4 mm prosthetic platform is connected to the implant and the grafting and suturing protocol from the Marginal Migration Concept follows. Following the same concept, the marginal portion of the first provisional

restoration will be supragingival in order to gain soft tissue in this critical area as we witness already 2 weeks post-op. 3 months post-surgery a second provisional restoration is fabricated in order to guide the maturation of the peri-implant soft tissue to the desired architecture. After 4 weeks second provisional restoration will be removed and an impression will be made in order to deliver the final restoration to the patient. The final restoration 2 weeks post cementation and 4 months post-surgery.

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Early Periimplant Tissue Healing on 1-Piece Implants With a Concave Transmucosal Design: A Histomorphometric Study in Dogs

Caroline Bolle, DDS, Marie-Paule Gustin, DS, PhD, Didier Fau, DVM, PhD, Patrick Exbrayat, DDS, Georges Boivin, PhD, and Brigitte Grosogoeat, DDS, PhD

Introduction: The purpose of our study was to investigate the early healing phase of marginal bone and soft tissues around unloaded 1-piece implants with a concave transmucosal design, in a dog model.

Methods: Twenty-four 1-piece implants with a concave transmucosal neck were inserted 1 mm subcrestally in the mandibular ridge of 8 beagle dogs. Four animals were sacrificed after 3 and 12 weeks of healing. Histomorphometric analysis was performed to measure the height of the periimplant tissues.

Results and Discussion: The overall height of the periimplant mucosa was, respectively, 2.67 and 2.52 mm, after 3 and 12 weeks. In the connective tissue, a soft tissue Oring seal was observed in the healing area provided by the transmucosal concavity, after 12 weeks. The location of the first bone-to-implant contact facing the implant shoulder was 0.00 and +0.18 mm, respectively, after 3 and 12 weeks of healing. Some bone apposition occurred on the implant shoulder during the healing.

Conclusion: Within the limits of the present study, a concave transmucosal design in 1-piece implants was associated with a short vertical value of biological width and promoted a mechanical interlocking of the implant body at the connective tissue and marginal bone levels.

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Rehabilitation of the Atrophic Posterior Mandible with Short (4mm) Implants: A case report

Roberto Pistilli, MD, Carlo Barausse, DDS, Luigi Checchi, MD, DDS, Pietro Felice, MD, DDS, PhD

This case report describes a successful implant-prosthetic rehabilitation of an atrophic posterior mandible with 4-mm-long implants. The patient refused to undergo any reconstructive surgery, and because the available bone up to the inferior alveolar nerve was only 5 mm or less, the patient received four implant of 4-mm length. Four months after implant placement, a provisional prosthesis was put in place; after another 4 months, this was then in turn replaced with a definitive prosthesis. The use of such short implants allows a fixed prosthetic solution without the need for vertically augmenting the mandibular bone. This procedure considerably reduces intra- and postoperative patient discomfort compared with reconstructive surgery for the placement of longer implants. The follow-up time was 1 year after implant loading.

Roberto Pistilli: Resident, Unit of Oral and Maxillofacial Surgery, A.C.O. San Filippo Neri, Rome, Italy.

Carlo Barausse: Resident, Dental School, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy.

Luigi Checchi: Full Professor, Dental school, unit of Periodontology and Implantology, University of Bologna, Bologna, Italy.

Pietro Felice: Assistant Professor, Dental School, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Int. J Periodontics Restorative Dent 2014; 34:713-718.
doi: 10.11607/prd.1733

Soft Tissue and Marginal Bone Adaptation on Platform-Switched Implants with a Morse Cone Connection: A Histomorphometric Study in Dogs.

Caroline Bolle, DDS, Marie-Paule Gustin, DS, PhD, Didier Fau, DVM, PhD, Patrick Exbrayat, DDS, Georges Boivin, PhD, and Brigitte Grosogeat, DDS, PhD

The purpose of this study was to investigate peri-implant tissue adaptation on platform-switched implants with Morse cone-type connection, after 3 and 12 weeks of healing in dogs. Ten weeks after mandibular premolar extraction, eight beagles dogs received three implants each. At each biopsy interval, four animals were sacrificed and biopsies were processed for histologic analysis. The height of the peri-implant mucosa was 2.32 mm and 2.88 mm, respectively, whereas the bone level in relation to the implant platform was -0.39 mm and -0.67 mm, respectively, after 3 and 12 weeks of healing. Within the limits of the present study, platform-switched implants exhibited reduced values of biologic width and marginal bone loss when compared with previous data..

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Int J Periodontics Restorative Dent 2016;36:221-228.
doi: 10.11607/prd.2254

Posterior jaws rehabilitated with partial prostheses supported by 4.0 x 4.0 mm or by longer implants: One-year post-loading results from a multicenter Randomised controlled trial.

Pietro Felice, MD, DDS, PhD, Luigi Checchi, MD, DDS, Carlo Barausse, DDS, Roberto Pistilli, MD, Gilberto Sammartino, MD, DDS, Irene Masi, DDS, Daniela Rita Ippolito, DDS

Purpose: To evaluate whether 4.0 x 4.0 mm dental implants could be an alternative to implants at least 8.5 mm long, which were placed in posterior jaw, in the presence of adequate bone volumes.

Materials and methods: One hundred and fifty patients with posterior (premolar and molar areas) jaws having at least 12.5 mm bone height above the mandibular canal or 11.5 mm below the maxillary sinus, were randomized according to a parallel group design, in order to receive one to three 4.0 mm-long implants or one to three implants which were at least 8.5 mm-long, at three centers. All implants had a diameter of 4.0 mm. Implants were loaded after 4 months with definitive screw-retained prostheses. Patients were followed up to 1-year post-loading and outcome measures were prosthesis and implant failures, any complications and peri-implant marginal bone level changes.

Results: Seventy-five patients were randomly allocated to each group. One patient dropped out after the 4-months post-loading evaluation from the long implant group. Up to 1-year post-loading, 3 patients lost one 4.0 mm-long implant each in comparison to 2 patients who lost one long implant each (difference in proportion = 0.013; 95% CI: -0.058 to 0.087; P=0.506). All failures occurred before loading, the failed implants were replaced and the delivery of two prostheses in each group was delayed for several months (difference in proportion = 0.0004; 95% CI: -0.068 to 0.069, P=0.685). Three short implant patients experienced three complications versus 2 long implant patients (difference in proportion = 0.013; 95% CI: -0.058 to 0.087; P= 0.506). There were no statistically significant differences in prosthesis failures, implant failures and complication. Patients with short implants lost on average 0.53 mm of peri-implant bone and patients with longer implants lost 0.57 mm. There were no statistically significant differences in bone level changes up to 1-year between short and long implants (mean difference = 0.038 mm; 95% CI: -0.068 to 0.138; P=0.198).

Conclusions: One year after loading 4.0 mm-long implants achieved similar results as 8.5 mm-long implants in posterior jaws, however 5- to 10-year post-loading data are necessary before reliable recommendations can be made.

Pietro Felice: Assistant Professor, Dental School, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Luigi Checchi: Full Professor, Dental school, unit of Periodontology and Implantology, University of Bologna, Bologna, Italy.

Carlo Barausse: Resident, Dental School, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy.

Roberto Pistilli: Resident, Unit of Oral and Maxillofacial Surgery, A.C.O. San Filippo Neri, Rome, Italy.

Gilberto Sammartino: Professor, Clinic of Oral surgery, University of Naples Federico II, Naples, Italy

Irene Masi: Resident, Department of Biomedical and Neuromotor Sciences, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy.

Daniela Rita Ippolito: Postgraduate Student, Department of Orthodontics, School of Dentistry, University of Brescia, Brescia, Italy.

Success rate of the EVL Evolution implants (SERF): A five-year longitudinal multicenter study.

Gérard Duminil, DDS, PhD; Michèle Muller-Bolla, DDS, PhD; Jean-Pierre Brun, DDS; Philippe Leclercq, DDS; Jean-Pierre Bernard, MD, PhD; David M. Dohan Ehrenfest, DDS, PhD

The purpose of this study was to evaluate the success rate of the SERF EVL evolution implants through a 5-year longitudinal multicentric study. Patients from 3 clinicians working in 3 different private practices (Grenoble, Nice, and Paris) and familiar with this implant system were included in this study; 413 patients and 1198 implants were followed over 5 years. The implant sites and implant types were recorded at the time of placement. The patients were followed yearly and controlled at the end of the study. The criterion for treatment evaluation or success was a qualitative variable related to 4 possible treatment outcomes: success, failure, ailing, and lost (dropout patients). Different variables (sex, bone quantity and quality at the implant site, location) were submitted to the chi-square test. A survival curve was established over 5 years according to the Kaplan Meyer method. The clinical follow-up was 3.1 6 1.2 years (ie, 1 to 6 years). At the end of this follow-up period, 1163 implants were classified as successful, 19 as failures, 12 as ailing, and 4 as lost (dropout). This implant system thus presented an overall success rate of 97.08%, over 5 years, independent of implant location, and for patient indications commonly encountered in private practice.

Gerard Duminil: DDS, PhD, is in private practice in Nice, France.

Michèle Muller-Bolla: DDS, PhD, is a professor, Department of Public Health, Faculty of Odontology at Nice University, Nice, France.

Jean-Pierre Brun: DDS, is in private practice in Grenoble, France.

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Jean-Pierre Bernard: MD, PhD, is a professor at the Department of Stomatology and Oral Surgery, School of Dental Medicine, University of Geneva, Geneva, Switzerland.

David M. Dohan Ehrenfest: DDS, PhD, is an assistant professor at Paris Descartes University, Paris.

Use of 3D cartilage scaffolds for the stabilization of implants and bone regeneration with fit-lock technique.

Giovanni Falisi¹, Massimo Galli¹, Pedro Vittorini-Velasquez², Juan C. Gallegos-Rivera³, Roberto Minasi¹, Alberto De Biase¹, Carlo Di Paolo¹

The surgical procedures for implant applications on the lateral upper areas depend on sinus pneumatization and availability of the residual bone. In these cases, autologous bone grafting remains the gold standard. Nevertheless, because of the morbidity associated to the donor site and the post-surgical complications, several alternative bone substitutes have been introduced, which, however, imply additional costs and show limited osteoinductive properties. Such limitations can be compensated with new regeneration strategies for biological and mechanical tissue restoration, a subject which has been addressed by tissue engineering in recent years. The authors present a new therapeutic option for implant application in the upper maxilla with bone availability less than 4 mm by using 3D scaffolds obtained from antigen-free porcine cartilage in the fit-lock technique. A longitudinal study on 18 consecutive cases was performed, with a 95.2% success rate one year after the implant. The advantages of this new technique are: 1) Functional and anatomical recovery of the maxillary antrum, 2) Immediate application of the implants; 3) Reduction of surgical times; 4) Absence of patient morbidity; 5) Local anesthesia; 6) Use of implants with a diameter > 4 mm.

Giovanni Falisi, Massimo Galli, Roberto Minasi, Alberto De Biase, Carlo Di Paolo: Department of Oral and Maxillo Facial Sciences, Sapienza University of Rome.

Pedro Vittorini-Velasquez: Maxillo-facial Clinic Cochabamba, Bolivia.

Juan C. Gallegos-Rivera: Department of Dentistry, Municipal Japanese University Hospital, Santa Cruz de la Sierra, Bolivia.

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
BIBLIOGRAPHIE

Becker J, Ferrari D, Mihatovic I, Sahn N, Schaer A, Schwarz F. Stability of crestal bone level at platform switched non-submerged titanium implants - a histomorphometrical study in dogs. J Clin Periodontol 2009; 36: 532-539.

Berglundh T, Abrahamsson I, Welander M, Lang NP, Linde H. Morphogenesis of the peri-implant mucosa : an experimental study in dogs. Clin Oral Impl Res. 18, 2007; 1-8

Hermann JS, Buser D, Schenk RK, Schoolfield JD, Cochran DL. Biologic width around one- and two-piece titanium implants. Clin Oral Implants Research 2001; 12: 559-571

Pearce AJ, Richards RG, Mitz S, Schneider E, Pearce RG. Animal models for implant biomaterial research in bone : a review. Eur Cell Mater 2007 Mar 2; 13: 1-10

Thomsen JS, Laib A, Koller B, Prohaskas S, Mosekilde LI, Gowin W. Stereological measures of trabecular bone structure: comparison of 3D micro computed tomography with 2D histological sections in human proximal tibial bone biopsies. Journal of Microscopy, vol 218, pt 2 May 2005, 171-179.

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Three concepts of single-unit restorations design to avoid biological complications

INTRODUCTION

The risk of biological complications has been rising with the number of implants placed every year. The restorative connection to the implant can be accomplished by either screwing or cementing the restoration on standardised or customised, individually shaped CAD/CAM abutments. Cement-retained prostheses are commonly used to restore single-unit restorations; however, there are reasonable disadvantages associated with a residual excess of cement after placing the crowns, which can lead to biological complications. Custom-made, screw-retained restorations enable precise modelling of gingival shape in accordance with the emergence profile, which provides optimal support for soft tissues and excludes the issue of undetected cement remnants.

PURPOSE

The aim of the study was to present three methods to restore a single implant unit that is activated by undetected residual cement rests while avoiding biological complications. Advantages of the presented concepts are discussed, based on the current literature and the utilization of the authors' own studies.

CLINICAL SITUATION

The concept No.1 is based on the individually designed and milled titanium abutment cemented extra-orally together with CAD/CAM full cerconia crown milled with the screw opening. The shoulder is on the level of the gingiva up to 1mm subgingivally for better esthetic outcome. The restoration is screw-retained.

Concept No.2 is a custom-produced abutment based on a standard titanium base cemented with a cerconia custom-milled part supporting the soft tissue and creating its margin 1 mm over the gingiva level. A full ceramic crown can be cemented on this type of custom abutment intra-orally for cement retained or cemented extra-orally with the screw opening for screw-retained restoration.

Restorative option No.3 is a full cerconia crown based on a standard abutment milled with the screw opening. The innovation of this crown compared to the standard procedure way of screw retained one is the custom-prepared subgingival part, which is marked on the gingiva level. The subgingival part was designed with respect to the emergence profile and made of unpolished cerconia due to perfect biocompatibility.

CONCLUSIONS

Decision making in restoring a single unit remains a challenge. Removal of all cement excess after crown cementation seems to be impossible, with the margin below the gingiva level, even using customized, anatomic abutments.

A few different approaches to manufacturing final restorations were introduced. The presented solutions combine advantages of screw-retained and cement-retained restorations. All of the reconstructions were manufactured according to the design of the emergence profile and biocompatibility, which can lead to predictable and safe conditions in the implant pocket over time.

Can implant design influence local tissue behaviour ?

PBI537

Thibaut MONEDIERE, Philippe DUCHATELARD, Jean-Luc VEYRUNE, Jean-François BOREL

Topic: Implant therapy outcomes, peri-implant biology aspects

The Twinkon® concept

- prosthetic crimping by its external morse taper (1)
- tissular crimping by its concave neck : the Tissue Creeping Profile (2)

Materials and methods


- A study conducted on 90 Twinkon® implants in 34 patients
- 5 years after their loading

Evaluation criteria :

- rate of peri-implant diseases : mucositis and peri-implantitis
- implant success according to Misch criteria
- radiological assessment of bone loss around these implants

Objectives


This study evaluates the long-term impact of the design of the Twinkon® implant on the stability of peri-implant tissues.

Radiological evaluation of marginal bone loss

Results

- 98,9 % survival rate
- 40,7 % mucositis
- 3,5 % peri-implantitis
- average bone loss near to zero (0,00 mm ± 0,72)
- bone stability or increased bone level despite some adverse clinical situations


Immediate implant placement

D0 D0 + 5 years

Thin biotype

D0 D0 + 5 years

Conclusion

The design of the Twinkon® implant must influence local tissue behaviour.

In certain situations, where local bone or mucous membrane-conditions didn't enable ideal three-dimensional implant positioning, we demonstrated peri-implant tissue-stability.

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PSA388

4 mm short implants in the treatment of reconstructive surgery failures of the atrophic posterior mandible: a retrospective study

R. Gasparro*, G. Sammartino*, I. Masi**, C. Barausse**, P. Felice**

* University of Naples Federico II, Department of Neuroscience, Reproductive Science and Odontostomatology, Naples, Italy
 ** University of Bologna, Department of Biomedical and NeuroMotor Sciences, Bologna, Italy

Topic: Implant therapy outcomes, surgical aspects

Background

Bone augmentation procedures can be associated with significant patient morbidity, can be more expensive and may require longer times (up to 1 year) before patients are rehabilitated with implant-supported prostheses. Moreover complications, especially for vertical augmentation techniques, are common.

Aim


The aim of this retrospective study was to evaluate whether short (4 mm) dental implants could be a suitable alternative in the failure of reconstructive surgery in the atrophic posterior mandible.

Materials and Methods

Patients experienced failure in the reconstructive surgery of the atrophic posterior mandible were enrolled. The type of reconstructive surgery, surgical complications, data related to the timing of failure and ways of resolution were collected. For their rehabilitation, super short (4 mm) implants (TwinKon® Universal System, Global D, Lyon, France) with a mean diameter of 4.09±0.19 mm were placed. After 4 months, provisional reinforced acrylic prostheses were delivered. Provisional prostheses were replaced, after 4 months, by definitive screw-retained metal-resin restorations. Outcome measures were prostheses and implant failures, any complications, peri-implant marginal bone level changes and patient satisfaction. Patients were followed up to an average of 2.37±0.98 years after loading.

Results

24 patients (5 Males and 19 Females; mean age: 60.5±7.4 years old) were enrolled in this study. They were treated with inlay technique, onlay technique, alveolar distraction osteogenesis and Guided Bone Regeneration (GBR). Complications were dehiscence/graft exposition, necrosis, peri-implantitis, paresthesia, non integrated graft and graft infection; they occurred on average 6.5±11.0 weeks after surgery. Significantly more complications occurred in the GBR group than the onlay and the distraction group (Kruskal Wallis test; P=0.001). 58 super short (4mm) implants were placed in the posterior mandible after the resolution of the complications (partial removal of the graft/membrane, total removal of the graft or irrigation with chlorhexidine followed by total removal of the graft). Up to a mean period of 2.37±0.98 years after loading, only two patients experienced short implant and/or prosthesis failures. The way of resolution of reconstructive complications did not influence the number of short implant complications. There were no statistically significant differences between marginal bone level changes at implant placement and at 1,2 and 3 years follow-up (P>0.05). The marginal bone loss at 3 years follow-up was 0.62±0.11 mm. All patients were fully satisfied with the implant treatment.

	Inlay (n=3)	Onlay (n=8)	Osteodistraction (n=4)	GBR (n=9)	Total (n=24)	P (chi-square test)
Dehiscence/graft exposition	3 (100%)	7 (87.5%)	0	6 (66.7%)	16 (66.7%)	0.013 *§
Necrosis	0	1 (12.5%)	4 (100%)	0	5 (20.8%)	0.000 *§²
Periimplantitis	2 (66.7%)	0	0	2 (22.2%)	4 (16.7%)	0.055
Non-integrated graft	0	1 (12.5%)	0	0	1 (4.2%)	0.625
Graft infection	0	0	0	4 (44.4%)	4 (16.7%)	0.067
Simultaneous graft exposition and infections	0	0	0	3 (33.3%)	3 (12.5%)	0.204

Table 1. Patient and intervention characteristics

Conclusions

This study suggests that super short (4 mm) dental implants could be a suitable, faster, cheaper and easier rehabilitation option in the failure of the reconstructive surgery of the atrophic posterior mandible. However, these results need to be confirmed by trials with follow-ups of at least 5 years.

References

- Esposito M, Pistilli R, Barausse C, Felice P. *European Journal of Oral Implantology*, 2014.
- Felice P, Cannizzaro G, Barausse C, Pistilli R, Esposito M. *European Journal of Oral Implantology*, 2014.

To access the posters online, please download the application

Dr. Pietro Felice DDS, MD, PhD
 Alma Mater Studiorum - University of Bologna
 Department of Biomedical and NeuroMotor Sciences
 pietro.felice@unibo.it

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