COM

Volume 2

PENDIUM

Clinique

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Partial maxillary rehabilitation with immediate provisional restoration - 8-year results


Mr. Rab, a 79-year old retired dental surgeon, visited our office because his cemented bridge at the maxillary has become mobile. The patient is healthy and there are no contraindications to implantology. The patient would like to have a fixed rehabilitation and to avoid having a removable prosthesis which he had a negative experience on the mandible.

We suggest the following treatment plan: removal of the remaining teeth excepted of his right maxillary bridge's. Placement of 6 Twinkon® implants according to the



immediate post-extraction technique and immediate loading with a screw retained provisional prosthesis, 24 hours after surgery.

6 months later, after tissues get mature, the final prosthesis is fabricated using scalloped abutments to preserve the peri-implant biologic width. At the abutment shift, we can observe the high quality of the mucosa thanks to the "tissue creeping profile" effect.


8 years later, clinical and X-ray examinations show a perfect bone and soft tissues stability.

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
Mandibular rehabilitation in CAD/CAM technique

A 66-year old Caucasian man came into our practice and was in good physical health with no medical risks. He had been wearing a complete denture on lower jaw and desired to rehabilitate his edentulous jaw with an implant-retained prosthesis.

A panoramic radiograph was taken and six In-Kone® implants were placed. We waited for a period of 3 months for the osseointegration process. In this case, due to the high aesthetic demand, the choice was made to fabricate a ceramic prosthesis.

In order to have the possibility to remove the rehabilitation, a screw-retained titanium implant bridge was proposed. A CAD/CAM solution provides the best precision and can be easily designed to obtain the best aesthetic outcome and to facilitate the oral hygiene manoeuvres in a more efficient and effective manner.

In conclusion, the framework is in titanium, the gingiva is in pink composite and the teeth have been processed in gold alloy-ceramic. The crowns were fabricated to be cemented on the framework, to facilitate the removal.

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
Complete upper and lower maxillary rehabilitation in immediate loading – 6-years results

A 54-year-old patient, a doctor and heavy smoker, comes to our office. He would like to restore his dentition with a fixed implants supported restoration and refuses to receive any removable alternative. The risk is important to comply with this demand even if the patient commits himself to improve the oral cares and reduce his cigarettes consumption. Yet, we accept.


The residual teeth are removed and 6 In-Kone® implants are placed respectively at the maxillary and mandible. Both arches are immediately loaded with two provisional bridges. Inter-prosthetic spaces are adapted to facilitate the maintenance.

4 months after surgery, the final prosthesis is fabricated, in ceramic at the upper jaw and resin at


the lower jaw. Considering the hygiene deficiency of this case, wider spaces should have been inserted, which the patient refused due to previous phonation problems.

6 years later, the oral hygiene conditions are still unfavorable but the X-rays don't show any inflammatory response. How to explain the resistance of this patient to peri-implantitis?

The taxons are certainly different from the periodontal disease's, the dento-gingival environment has been modified, the teeth have been replaced by a smooth and sealed implant connexion and the auto-immune function of the patient is certainly performing well.

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

Twinkon® single unit restoration – 3-month and 7-year results

A 70-year old patient comes to our office to replace his missing tooth on 24. The latter has been removed two months ago after it fractured. The case is treated with a Twinkon® implant: a transucosal implant with a concave machined collar which purpose is to create a strong mucosal joint and to stabilize the alveolar bone in the long term.

Let us note that the patient already received a Twinkon® implant 7 year ago (former design with a single thread) to restore a single crown on 26. 7 years later, the results on this patient seems to comfort the relevance of this concave profile with regards to the bone-mucosa complex stability.

The patient to be treated presents an important class II overbite. The aesthetic demand of this second restoration in position 24 requires a precise tridimensional positioning of the future emergence profile. The drilling protocol we use is the ULTIMATE one which corresponds to the new Twinkon® design (twinthread) and which allows a uniform distribution of the primary stability over the bone-to-implant contact surface.

A single unit ceramic crown is fabricated and cemented. It will be interesting in the long term to compare the surrounding bone behavior between the two implant bodies designs.

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Complex multiple edentulous restoration - 5-year results

The edentulated patient consulted in April 2010. The difficulty of this case can be appreciated on the clinical view (picture 1) as well as the X-ray (picture 2). The study cast is mounted in a semi-adjustable articulator to define the treatment plan and to ease the communication with the patient (picture 3). On the maxillary, the treatment includes a bridge from 13 to 23, two crowns on 24 and 25, one implant on 26 and an implant-supported bridge from 14 to 16. On the mandible, it includes a bridge from 46 to 33 on natural teeth and an implant-supported bridge from 34 to 36.

First, a bilateral sinus lift is performed in August 2010. Provisional bridges are fabricated from the articulated study models (picture 4). The bridges are placed after the natural teeth have been prepared on both arches (pictures 5-6-7). Five implants are placed in May 2011 (picture 8). The final prostheses are placed in September 2011 (photos 9-10). The situation on October 2016, 5 years later (pictures 11-12)

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

Maxillary complete denture rehabilitation with guided surgery and immediate loading


A 50 year-old edentulous patient, ASA 1, visited our office to have a fixed prosthetic restoration with minimal post-operative sensitivity. We propose an "All on Four" type rehabilitation together with a guided surgery protocol in order to reduce the chair time and the wounds due to the flap thanks to a minimally invasive flapless surgery technique

The stereolithic surgical guide allows us to manage an optimal positioning of the implant with regards to the anatomy of the patient, the bone volume and the emergence of the final prosthesis abutments. The main steps of our protocol are:

- ConeBeam CT with an implant carrier equipped with radiopaque Pixmarkers
- Transmission of the DICOM data and validation of the "All on Four" planning

- Fabrication of a surgical guide, an occlusal guide and a custom open impression tray
- Placement of 4 implants via surgical guide (flapless technique)
- Selection and tightening of the conical abutments
- Registration of intermaxillary positions with the occlusal guide
- Impression taking with the custom open impression tray
- Tightening of the screw retained bridge.

In conclusion, the pre-planning allows to optimize the 3D placement of the implants with regards to the prosthetic project and the bone volume while providing security to a minimally invasive surgical approach.

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
Single unit premolar restoration and 12-month follow up

A 63-year old non-smoker patient comes to our office further to an apical lesion on 14. The tooth is removed, and the site sutured.

3 months later, a Ø3.5 L11.5 mm In-Kone® implant is placed in one surgical time.


After a three months period again, we proceed

to the fabrication of the zirconia screw retained prosthesis.


The X-ray control at 12 months shows the good health of the alveolar bone as well as the peri-implant mucosa.

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
Reconstruction of an alveolar ridge with allogenic bone to place implants - 2-year results

A 35 year old female patient comes to our office for she suffers an aesthetic failure further to a misplaced implant on 22 as well as several infection episodes.

An X-ray examination shows a large cystic zone spreading from the root of the 21 to the mesial side of the implant. Another cyst appears also at the tip of the 11.

The proposed treatment plan is the following: the implant must be removed, as well as 21 and 11 whose remaining value is not sufficient.

A temporary removal prosthetic is placed after the surgery. After two months of healing, a bone grafting is performed to rebuild the ridge.

The patient refuses that her chin bone nor her

retro molar bone be harvested even though they are conventional protocols. So, we decide to use a cortico-spongyous allogenic bone block graft from BIOBank®. It is placed on the buccal side of the alveolar ridge with a traditional surgery procedure. After four months of healing, three In-Kone® implants are placed into a suitable bone volume.


Temporary screw-retained crowns are adjusted on the implants three months and a half later to manage the surrounding soft tissues and prepare the emergence profile prior to the impression taking and final prosthetics fabrication.

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



CAD/CAM implant bar for fixed-removable restoration

A 63-year old patient visits our clinic because he complains about his current prosthesis' mobility which affects his chewing ability. The maxillary anterior teeth are loose (severe loss of attachment over 80%). The teeth must be removed and replaced by a maxillary fixed-removable prosthesis.


4 In-Kone® implants Ø4 L11,5 are placed to support an immediate provisional prosthesis the day of the surgery. Two months later, a second surgery is performed to place conical abutments. A pick-up impression is taken with an open tray. The occlusal registration is achieved with a duplicate of the immediate prosthesis which is screw-retained onto the implants (in position 13 and 23) with fixed castable sleeves.

A master model with gingival mask is casted by the technician in order to fabricate a CAD/CAM bar. This bar provides an optimal distribution of the Locator® and forces onto the 4 implants. The technician fabricates a metal base plate which can accept Locator® removable attachment.


The patient can thus enjoy both an aesthetic and functional prosthesis. The bar, the abutments emergence and the overdenture intrados are easy to clean. The benefits of this protocol are the good distribution of the forces and the maintenance accessibility which both contributes to the long-term success of the rehabilitation.

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Complete maxillary and mandibular rehabilitation with immediate loading


This patient comes to our office and introduces himself as a dynamic and busy entrepreneur, carrier minded and running short of time to take care of his personal health issues. He would like a complete restoration of his smile and masticatory ability.

The initial X-ray shows numerous non-treated teeth avulsions in an advanced destructive chronic periodontitis environment. After a preliminary treatment period of the disease of two weeks, we evaluate again the situation. The patient has been involved in the treatment and positively reacted. However, the restoration of the remaining teeth is likely to be a short term solution which cannot meet the patient initial requirements.


Thus, we decide to restore the patient under sedation with a complete maxillary and mandibular implant supported in one step. After the remaining teeth are removed, 10 and 8

implants are placed respectively at the maxillary and mandible, without surgical guide. The bone defects are grafted with vestibular allogenic bone graft from BIOBank® + A-PRF and a “Steaky Bone” appositional graft + I-PRF on both arches. Partial thickness flaps going beyond the grafting area are coronally repositioned and sewed up without tension in “apical mattress” and single papilla sling suture.

A 48 hour post-operative provisional bridge is placed at the lower jaw. 48 hours later, the upper jaw is restored with a similar provisional bridge. 6 months later, the final restoration (Zircon – Ceramic) is fabricated and placed. It is interesting to observe the quality of the keratinized gingiva as well as the presence of bone supported interdental papilla. This clinical result can be obtained thanks to sub-crestal “platform-switched” implants combined with bacteria sealed morse taper prosthetic interfaces.

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Partially posterior rehabilitation 2-year results

The patient is a 48-year old woman, who visit our clinic and complains about discomfort due to the partially posterior edentulous space of her mandible. The patient would like to restore this defect with a prosthesis providing maximum biomimetic properties.

Clinical and X-ray examinations show a vertical bone loss which makes impossible to place implants likely to meet her demand without going through bone augmentation. Accordingly, the patient is informed about various surgical techniques options of bone augmentation in order to create the suitable volume for implants placement. The patient decides to get a bone augmentation using a customized allogenic bone block technique.

Based on patient's mandible CT scan, pre-

planned allogenic bone block is customized with CAD/CAM technology. The bone block is fixed to the alveolar ridge with Graftek[®] screws and covered with a thin layer of xenogeneic bone graft granules plus a resorbable collagen membrane. 6 months later, the augmented site re-entry is made and three In-Kone[®] UNIVERSAL dental implants are placed in 44,45,46, 4 months later, customized zirconium abutments on titanium bases and ceramic dental restorations are fabricated and placed on the implants.

The result of the dental treatment fully meets the expectations of the patient. Clinical and X-ray observations one and two years later show healthy peri-implant soft tissue in the augmented site, integrity of implant-based restorations, a stable functional and aesthetic result.

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







Rehabilitation of the atrophic posterior mandible with 4 mm short implants

A 46-year-old patient visited our clinic for a treatment to replace the missing teeth 45, 46, 47. Clinical and X-ray examinations showed a significant resorption of the alveolar bone ridge. The amount of bone available over the mandibular nerve does not exceed 5 mm. Accordingly, the patient was proposed several dental treatment options including bone augmentation prior to long implants placement, 4 mm short implants placement or removable prosthesis. The patient compromised with the short implants treatment option because he fears invasive surgery and doesn't accept to have a removable prosthesis.

Short implants Twinkon® 4 have been chosen based on their successful clinical application and scientific publications that demonstrate their predictability. Three Twinkon® 4 implants (Ø4.5 L4 mm) have been placed in 45, 46,

47 under local anesthesia. The Twinkon® 4 implants have been carefully positioned with regards to the prosthetic project, in a single-stage surgery. Healing screws were removed from the implants 4 months after surgery and conical abutments have been installed. A screw-retained bridge has been fabricated from zirconium dioxide ceramic based on titanium abutments. An X-ray control confirms the osteointegration of the implants as well as the fitting of the prosthetic interfaces.

Considering the significant bone defect, the prosthesis treatment has been feasible in only for 4 months thanks to short implants. Functional and aesthetic results meet both the clinician and patient expectations.

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

Rehabilitation of a double maxillary agenesis

A 52-year-old ASA1 female patient suffers a double agenesis in 12 and 22 which was restored with a Maryland bridge. She doesn't like her smile.

The tridimensional X-ray and scanner evaluations of the edentulous sites reveals an available bone crest of 2 to 3 mm as well as converging apex of the neighboring canines and premolars.

The treatment we propose includes the preparation of the interdental spaces with ali-


gners both at the apical and coronal levels.

Then, a double guided split is performed and followed by the placement of two narrow implants 3.0 mm.


After tissue maturation, the e-Smile laboratory fabricates in CAD/CAM technique full ceramic crowns to be cemented on activated angulated abutments.

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





Inferior alveolar nerve transpositioning for dental implant rehabilitation - 4-year results

A 59-year-old male presented jaw bone defect caused by gunshot that iliac block bone graft was done 40 years ago. His chief complaint was unable to eat with the traditional denture he was using over 40 years.

Intraoral examination revealed serious mandibular bone loss so that the alveolar crest was in the same level with the floor of the mouth. No keratinized tissue was observed. Radiographic diagnosis showed the inferior nerve crossed the midline and was close to the alveolar crest. It even exposed on the crest in the position 41, 31. The residual bone height was 7 to 11 mm. In such a situation, it's impossible to perform block bone graft, bone distraction or All on 4.


Our decision was transpositioning inferior alveolar nerve in combination with vertical iliac

block graft and 4 In-Kone® implants placement simultaneously. Implant obtained good initial stability ($\geq 35\text{N/cm}$). 2 months later, we carried out free soft tissue graft and vestibuloplasty to improve soft tissue around the implants and loose connective tissue on alveolar crest.


A bar-retained implant overdenture with clip attachment was delivered. However, after 4 years loading, the overdenture was unstable and unretentive as the clip worn out. The new bar-retained overdenture with equator attachment from Rhein'83 system was replaced. Good osseointegration around implants are still maintained with no bone loss after 4 years. The patient is satisfied with the new overdenture that improve significantly his chewing function and aesthetic, especially his quality of life.

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





Rehabilitation of a lateral incisor

A 47-year old female patient consulted because she is worried about the mobility of her lateral incisor (12). The X-ray shows a root fracture. After the extraction and curettage of the socket, a buccal cortical bone fenestration is observed. A bone and connective tissue graft is operated and protected by a temporary Maryland bridge. 6 months later, an In-Kone® implant


is placed and loaded with an immediate nonfunctional provisionalization. After the soft tissue maturation phase, an Emax screw retained single restoration is fabricated (LDL technician).

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
4 mm-long versus longer implants in augmented bone in posterior atrophic jaws: 4-month post-loading results from a multicenter randomised controlled trial

Marco Esposito, Giovanni Zucchelli, Carlo Barausse, Roberto Pistilli, Anna Trullenque-Eriksson, Pietro Felice

Purpose: To evaluate whether 4-mm long dental implants could be an alternative to augmentation with equine bone blocks and the placement of at least 10-mm long implants in atrophic posterior jaws.

Materials and methods: Forty patients with atrophic posterior (premolar and molar areas) mandibles having 5 to 6 mm bone height above the mandibular canal and 40 patients with atrophic maxillae having 4 to 5 mm below the maxillary sinus, were randomised according to a parallel group design to receive one to three 4.0 mm-long implants or one to three implants, which were at least 10 mm long, in augmented bone at two centres. All implants had a diameter of 4.0 or 4.5 mm. Mandibles were vertically augmented with interpositional equine bone blocks and resorbable barriers. Implants were placed 4 months after interpositional grafting. Maxillary sinuses were augmented with particulated porcine bone via a lateral window covered with resorbable barriers, and implants were placed simultaneously. Implants were not submerged and were loaded after 4 months with provisional prostheses. Four months later, screw-retained reinforced acrylic restorations were delivered, and then replaced after 4 months by definitive screw-retained metal-composite prostheses. Patients were followed up to 4-months post-loading. Outcome measures included prosthesis and implant failures, any complication and peri-implant marginal bone level changes.

Results: No patient dropped out. In six augmented mandibles (30%), it was not possible to place implants which were at least 10.0-mm long, therefore shorter implants had to be placed instead. In particular, one mandible fractured and the patient did not want to go ahead with the treatment. One implant of the patient with the mandible fracture from the augmented group failed versus two 4.0 mm implants in two patients from the short implant group. In the maxillae, three short implants failed in two patients versus five long implants in three patients (two long implants and one short implant dropped into the maxillary sinus). Two prostheses on short implants (one mandibular and one maxillary) were placed at a later stage because of implant failures versus four prostheses (one mandibular and three maxillary) at augmented sites. In particular, three patients of the augmented group (one mandible and two maxillary) were not prosthetically rehabilitated. There were no statistically significant differences in implant failures (P (chi-square test) = 1.000; difference in proportion = 0; 95% CI: -0.13 to 0.13) or prostheses failures (P (chi-square test) = 0.399; difference in proportion = 0.05; 95% CI: -0.06 to 0.16). At mandibular sites, nine augmented patients were affected by complications versus one patient treated with short implants (P (chi-square

test) = 0.003; difference in proportion = 0.40; 95% CI: 0.16 to 0.64), with the difference being statistically significant. No significant differences were found for the maxillae: eight sinus lift patients versus three patients rehabilitated with maxillary short implants were affected by complications (P (chi-square test) = 0.077; difference in proportion = 0.25; 95% CI: -0.02 to 0.52). Patients with mandibular short implants lost on average 0.40 mm of peri-implant bone at 4 months and patients with 10 mm or longer mandibular implants lost 0.52 mm. Patients with short maxillary implants lost on average 0.48 mm peri-implant bone at 4 months and patients with 10 mm or longer maxillary implants lost 0.50 mm. The difference was statistically significant in the mandibles (mean difference: -0.12 mm, 95% CI: -0.20 to -0.04, P (ANCOVA) = 0.006), but not in the maxillae (mean difference: -0.02 mm, 95% CI: -0.10 to 0.07, P (ANCOVA) = 0.711).

Conclusion: Four months after loading 4.0 mm-long implants achieved similar results, if not better, than longer implants in augmented jaws, but were affected by fewer complications. Short implants might be a preferable choice to bone augmentation, especially in mandibles, since the treatment is less invasive, faster, cheaper, and associated with less morbidity; however, 5- to 10-year post-loading data is necessary before making reliable recommendations..

Marco Esposito DDS, PhD: Freelance researcher and Associate Professor, Department of Biomaterials, The Sahlgrenska Academy at Göteborg University, Sweden

Giovanni Zucchelli DDS, PhD: Associate professor, Department of Biomedical and Neuromotor Sciences, University of Bologna, Bologna, Italy

Carlo Barausse DDS: Resident, Department of Biomedical and Neuromotor Sciences, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Roberto Pistilli MD: Resident, Oral and Maxillofacial Unit, San Camillo Hospital, Rome, Italy; Private practice, Rome, Italy

Anna Trullenque-Eriksson DDS, PhD: Folkandvården Sylte, Trollhättan, Sweden; Honorary Research Fellow, Institute of Dentistry, Barts and The London School of Medicine and Dentistry, Queen Mary University of London, London, UK

Pietro Felice MD, DDS, PhD: Researcher, Department of Biomedical and Neuromotor Sciences, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Implant stability evaluation by resonance frequency analysis in the fit lock technique. A clinical study

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Summary: Surgical procedures for the application of implants in the lateral-superior sector or affected availability of the residual bone.

When this condition is lower than 5 mm it is recommended, the technique involving two therapeutic phases, a reconstructive and an applicative one, as reported in the international literature, are adopted. The authors propose here a new method with the potential to apply implant simultaneously with the reconstructive phase.

The aim of this longitudinal retrospective study was to evaluate the stability of implants applied with the fit lock technique in the upper maxillary's sinus with bone availability lower than 4 mm by measuring resonance frequency at different follow-up periods. The measurement, carried out on 30 implants, were analyzed with specific statistical procedures.

The results indicate the stability of the implants inserted with the fit lock method increases progressively over time in a statistical significant manner. The stability recorded after one year from the insertion (ISQ T2) is significantly higher than the recorded after six months (ISQ T1), and this is significantly higher than that recorded at the time on implants placement (ISQ T0)

The implants inserted in the maxillary zones with scarce bone availability and applied with this technique showed a similar stability as reported with other techniques. In light of the results, the authors confirm that the primary stability represents the basic requirements to guarantee a correct healing of the implants and demonstrates that the fit lock technique also owes reaching these conditions when bone availability is minimal

Key words: resonance frequency analysis, bone grafting, dental implants.

Giovanni Falisi, DDS, PhD: Department of Oral and Maxillo-Facial Sciences University of Rome, Italy

Massimo Galli, MDS: Department of Oral and Maxillo-Facial Sciences University of Rome, Italy

Perdo Vittorini Velasquez, MDS: Maxillo-Facial clinic Cochabamba, Bolivia

Juan Carlos Gallegos Rivera, DDS: Department of Dentistry, Municipal Japanese University Hospital, Santa Cruz de la Sierra, Bolivia

Carlo Di Paolo, MDS: Department of Oral and Maxillo-Facial Sciences University of Rome, Italy

Early loading versus immediate loading: case examples

Alberto Maltagliati, Andrea Ottonello, Giulio Raffaghello, and Andrea Mascolo

Introduction: “Immediate loading” is a favorite way to maintain appropriate anatomic and topographic proportions between dental arches, to maintain the contact with periodontal tissues, and for the “occlusal / contact memory”.¹⁻² An edentulous area that remains this way for a long period of time may modify and even preclude the possibility of recreating an occlusion that is esthetic, functional, always predictable, and similar to the original.³ The concept of “primary stability,” torque wrench, osseointegration, and adequate/suitable denture, are the fundamental factors for immediate and long-distance clinical results.⁴⁻⁵ With immediate loading, the management of peri-implant soft tissues can lead to a long-term predictability with 30-year old case-studies;⁶⁻⁷ moreover, the waiting time for osseointegration allows for modulation/mediation of some human mistakes not managed in the early loading phases in the 72-hour prosthetic technique.⁸ The use of autologous fibrin glue in the same surgery increases the neoangiogenesis and the healing factors’ activity, also obtaining a better predictability regarding soft-tissue management.⁹⁻¹⁰

Aim of the study: The aim of this study was to demonstrate the concept of middle and long-term esthetic and functional predictability differentiating: - full arch implant restoration, cemented technique, with loading after 60-90 days (early loading) - full arch implant restoration, screwed technique, with loading after 72 hours (immediate loading) since fixture insertion.

Conclusions: Finally, in our opinion, supported from the recorded data, the immediate loading technique (full-arch with screwed technique) represents the current “gold standard” in implant prosthetic restoration, because it reduces intraoperative time and increases wound healing, in hard and soft tissues, with excellent patient compliance.

Alberto Maltagliati, DDS: Degree from the University of Genoa, in Italy. He has taught at the University of Genoa and the University of Chieti.

Andrea Ottonello, DDS, PhD: Degree from the University of Genoa, and also earned a PhD at the same university. PhD student and Fellowship in Periodontology at the Department of Periodontics and Oral Medicine at the University of Michigan. Professor in Periodontology at the University of Genoa.

Giulio Raffaghello, DDS: graduated from the University of Pavia, Italy, and attended a course in endodontics. Private practice, with Dr. Mascolo, DDS, in research activities.

Andrea Mascolo, DDS: has a Master in Oral Surgery from the University of Pisa, Italy. He qualified for Implantology at the Brånemark Clinic, Gothenburg, Sweden and attended the Summer School Clinical Implantology Periodontology at the University of Heidelberg, Germany. Founder and active member of the Computer Aided Implantology (CAI) Academy. Professor at the University of Genoa since 2001. Speaker and author in Italy and abroad, and is an international peer reviewer. Private practice in Novi Ligure, Italy, devoted to minimally invasive surgical procedures on compromised patients. UK GDC registered, and is an expert in periodontics and implantology in private clinics.

The effects of surgical preparation techniques and implant macro-geometry on primary stability: An in vitro study

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Abstract

Background: The attainment of a good primary stability is a necessary condition to ensure the success of osseointegration in implantology. In type IV cancellous bone, however, it is possible that a reduced primary stability can lead to an increased rate of failure. The aim of this study was therefore to determine, with the help of the resonance frequency (Osstell mentor), which technique of implant site preparation (piezo surgery, conventional, under-preparation, bone compaction, osteodistraction) and macro-geometry is able to improve implant stability in type IV cancellous bone.

Material and Methods: 10 pig ribs were prepared with a surgical pre-drilled guide, calibrated for a correct implant positioning. On each rib, 5 implant sites (one for each technique) were prepared. Successively, 50 conical implants (Tekka Global D) were inserted and measured with the resonance frequency to evaluate the primary stability. Data collected were analyzed by analysis of variance (ANOVA) to test whether the Implant Stability Quotient (ISQ) values of the five techniques were significantly different.

Results: The results showed that no significant differences among the ISQ values of the five techniques used were found. Also, no significant differences in the macro-geometry of the two types of compared implants were observed. However, the macro-geometry of Tekka implants, characterized by a double condensing thread, seems to provide greater ISQ values than those of single thread implants when using the same technique.

Conclusions: In light of these preliminary data, it is conceivable that in cases of reduced stability, such as those occurring with a type IV bone, all means ameliorating the primary stability and accelerating the osseointegration can be utilized.

Keywords: Implant primary stability, Resonance frequency analysis, Implant site preparation

Giovanni Falisi, MD: Department of Life, Health, and Environmental Sciences, School of Dentistry, University of L'Aquila, L'Aquila, Italy
Marco Severino, MD: Department of Life, Health, and Environmental Sciences, School of Dentistry, University of L'Aquila, L'Aquila, Italy
Claudio Rastelli, MD: Department of Life, Health, and Environmental Sciences, School of Dentistry, University of L'Aquila, L'Aquila, Italy
Sara Bernardi, PhD: Department of Life, Health, and Environmental Sciences, School of Dentistry, University of L'Aquila, L'Aquila, Italy
Silvia Caruso, DDS: Department of Life, Health, and Environmental Sciences, School of Dentistry, University of L'Aquila, L'Aquila, Italy
Massimo Galli, MD: Department of Oral and Maxillofacial Sciences, "Sapienza" University of Rome, Italy
Luca Lamazza, MD: Department of Oral and Maxillofacial Sciences, "Sapienza" University of Rome, Italy
Carlo Di Paolo, MD: Department of Oral and Maxillofacial Sciences, "Sapienza" University of Rome, Italy

Four-year post-loading results of full-arch rehabilitation with immediate placement and immediate loading implants: A retrospective controlled study

Simonpieri, Alain / Gasparro, Roberta / Pantaleo, Giuseppe / Mignogna, Jolanda / Riccitiello, Francesco / Sammartino, Gilberto

Objectives: The aim of this study was to evaluate implant survival rate and to measure peri-implant bone changes in full-arch rehabilitations with immediate placement and immediate loading implants with platform switching and Morse taper connection, in addition to platelet-rich fibrin (PRF) and buccal bone augmentation, after 4 years of follow-up.

Method and Materials: In this retrospective controlled study, patients who had been fully rehabilitated with immediate placement and immediate loading implants were evaluated 4 years post-loading. Implants with platform switching and Morse taper connections were used (In-Kone Universal System, Global D) and PRF and buccal bone augmentation were applied. The radiographic bone loss was calculated by subtracting the bone level at baseline (BLT0) from that at the 4-year follow-up (BLT4) in immediate and delayed implants. Measurements were made at the distal, mesial, vestibular, and oral sites of the implants and the deepest value was recorded. Implants placed in extraction sites and implants placed in healed sites were considered. A comparison between the groups was performed using the Mann-Whitney test. The implant survival rate was calculated using the Kaplan-Meier analysis.

Results: In total, 42 patients (28 females and 14 males; average age 55.8 years old, age range 45 to 77) were recruited to this study. A total of 334 implants were put in place (226 in the maxilla, 108 in the mandible). The implant survival rate was 97.8% for the maxilla and 98.1% for the mandible, 98.3% for immediate implants and 96.9% for delayed implants. No statistically significant differences ($P > .05$) in the mean radiographic bone loss (mBL) were observed when comparing the immediate and delayed implants and the anterior and posterior implants. Statistically significant differences were found in the mBL between the mandibular and maxilla implants in the vestibular ($P = .01$) and mesial ($P = .001$) sites.

Conclusion: Within their limits, the present results suggest that rehabilitation with platform switching and taper connection implants, in addition to buccal bone augmentation and the use of PRF, can lead to predictable results.

Keywords: bone loss, dental implant, dental implant-abutment connection, immediate dental implant loading, survival rate

A Minimally Invasive Approach Using a 4-mm Implant without extraction of Impacted Maxillary Canine: Four-Year Postloading Results

Pietro Felice, MD, DDS, PhD, Carlo Barausse, DDS, Martina Stefanini, DDS, PhD, Roberto Pistilli, MD, Giovanni Zucchelli, DDS, PhD

The aim of this case report was to suggest an alternative minimally invasive surgical approach to an impacted maxillary canine using a 4-mm-long implant for a fixed prosthetic rehabilitation, avoiding tooth extraction or surgically forced extrusion and exploiting the 6 mm coronal bone availability. At 4 years postloading, the implant was healthy and well integrated with stable marginal bone levels. The 4-mm length of the implant reduced operative times, postsurgical morbidity, possible complications, and costs. Short implants might be an alternative to traditional, more invasive surgical procedures used in the rehabilitative treatment of impacted maxillary canines.

Pietro Felice MD, DDS, PhD: Researcher, Department of Biomedical and Neuromotor Sciences, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Carlo Barausse DDS: Resident, Department of Biomedical and Neuromotor Sciences, Unit of Periodontology and Implantology, University of Bologna, Bologna, Italy

Martina Stefanini, DDS, PhD: Resident, Department of Biomedical and Neuromotor Sciences, University of Bologna, Bologna, Italy

Roberto Pistilli, MD: Resident, Oral and Maxillofacial Unit, San Camillo Hospital, Rome, Italy

Giovanni Zucchelli DDS, PhD: Associate professor, Department of Biomedical and Neuromotor Sciences, University of Bologna, Bologna, Italy

Publication pending in Int J Periodontics Restorative Dent – dec. 2017

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